


# Spenco Kids® PolySorb® INSOLES

Replacement Insoles



Cat. No. 45-039

## Memory foam comfort and arch support.

### Target Consumer:

Consumers with tired feet.  
Consumers needing additional cushioning and shock absorption for feet stress.

### Recommended For:

Reducing pain associated with Sever's heel, heel spurs, and foot impact.  
Superior cushioning and comfort.  
Energy return helps reduce fatigue for better performance.



### Best For Use In:

Athletic, tennis, basketball and training shoes with removable insoles.



### Components\*:

Memory Foam, EVA, SpenCore® Material, PU Gel, Fabric, Thickness at Heel: 12.50mm (1/2in), Thickness at Forefoot: 6.00mm (1/4in)



### Instructions:

Remove the original insole from the shoe (some brands of shoes glue the insole to the shoe, but careful removal will not damage the shoe). If necessary, use the original insole as a guide to trace and trim. Insert the Spenco® Insole into the shoe.



### Care:

Hand wash in warm water. Air dry.

### Sizing Chart:

Size	0	1	2
Kid's	3/4 1/2	5/6 1/2	7/8

### Features & Benefits:

- 1 **Memory Foam** Adjusts for customized cushioning for the entire foot.
- 2 **Molded Polyurethane Base** absorbs shock and supports the arch and heel.
- 3 **4-way Stretch Fabric with Silpure Antimicrobial** helps prevent blisters while controlling odor.
- 4 **SpenCore® Material** provides the ultimate in long lasting cushioning through spring-like shock absorption.
- 5 **Additional Heel Strike** provides maximum cushioning at heel strike.

### One-Year Unconditional Guarantee

Call toll-free: 1-800-877-3626  
In Canada: 1-800-387-9538

 [www.spenco.com](http://www.spenco.com)



\*Measurements are approximate and vary by size.  
Silpure is a Registered Trademark of Thomson Research Associates, Inc.