


# Spenco® GEL HEELCUPS

Gel Heel Cups



Cat. No. 39-826

## Deep cupping and gel cushioning.

### 🎯 Target Consumer:

Consumers seeking contoured cushioning around the heel area.

Consumers suffering from general heel pain, heel spurs or bruises.

### ✅ Recommended For:

Daily wear – absorbing shock to reduce heel pain.



### Best For Use In:

Any casual, work or walking shoe.



### Components\*:

TPR Gel (dual density),

Heel thickness – 9.22mm [3/8"]



### Instructions:

Place Spenco® Gel Heel Cup against the inside heel area of shoe. For best fit, unlace shoe and allow maximum room to correctly position foot in the heel cup. Tighten laces after foot is comfortably in place.



### Care:

Hand wash with warm water. Air dry.

## Features & Benefits:

1

**TPR Comfort Gel** Thermoplastic Rubber combines superior cushioning and energy return for ultimate comfort.

2

### Anatomical Cushioning System

Dual density gels absorb impact and shock where the heel needs it most.

3

**Heel Cup** provides extra cushioning for the heel.

## One-Year Unconditional Guarantee

COMFORT

## Sizing Chart:

Size	Small/Medium	Medium/Large
U.S. Women's	5-8	9-12
U.S. Men's	6-7	8-15
U.K.	4-6	7-14
European	36-40	41-48

\*Measurements are approximate and vary by size.

Call toll-free: 1-800-877-3626  
In Canada: 1-800-387-9538

 [www.spenco.com](http://www.spenco.com)

