


Spenco® GEL

TOTALSUPPORT

Total Support Gel Insoles



Cat. No. 46-300

Total support and gel cushioning.

Target Consumer:

Consumers seeking all day comfort for their feet.
Consumers standing for long periods of time:
workers, students and travelers.

Recommended For:

Daily wear – absorbing shock to reduce foot fatigue
and sore feet.
For workers, acts like an anti-fatigue mat in the shoe.



Best For Use In:

Any casual, work or walking shoe.



Components*:

Heel thickness – 8.00mm [5/16"],
Forefoot thickness – 5.00mm [3/16"],
Thickness at toe – 3.00mm [1/8"],
Material – TPR Gel,
Fabric – Polyester w/Silpure Antimicrobial



Instructions:

Remove existing insole (some brands of shoes glue
the insole to the shoe, but careful removal will not
damage the shoe). Insert the Spenco® Gel Insole into
the shoe. Trim with scissors if necessary.



Care:

Hand wash with warm water. Air dry.

Features & Benefits:

- 1 **TPR Comfort Gel** thermoplastic rubber combines superior cushioning and energy return for ultimate comfort.
- 2 **Metatarsal Arch Support** relieves pressure from the ball-of-foot.
- 3 **Top Cloth with Silpure Antimicrobial** helps eliminate odor causing bacteria and reduces friction.
- 4 **Stability Cradle** supports the arch and heel for greater comfort. Makes it easier to place into the shoe.
- 5 **3-POD™ Cushioning System** absorbs impact at heel strike.

One-Year Unconditional Guarantee

Sizing Chart:

Size	0	1	2	3	4	5	6
U.S. Women's	3/4	5/6	7/8	9/10	11/12		
U.S. Men's			6/7	8/9	10/11	12/13	14/15
U.K.	2/3	4/5	5/6	7/8	9/10	11/12	13/14
European	32/34	36/38	38/40	40/42	42/44	44/46	46/48

Call toll-free: 1-800-877-3626
In Canada: 1-800-387-9538

 www.spenco.com



*Measurements are approximate and vary by size.
Silpure is a Registered Trademark of Thomson Research Associates, Inc.